

The natural way of Sweden

# SISSEL® ACUPRESSUR MAT

- | Gebrauchsanweisung
- | Gebruiksaanwijzing
- | Mode d'emploi
- | Instructions for use
- | Instrucciones de uso
- | Istruzioni per l'uso
- | Bruksanvisning
- | Käyttöohje

# SISSEL<sup>®</sup> ACUPRESSUR MAT

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

# SISSEL® ACUPRESSUR MAT

## BLISSFUL RELIEF FOR TENSE AND

Our Swedish SISSEL® Acupressur Mat is based on ancient Indian traditions. It has a blissful and relaxing effect in the treatment of neck, shoulder and back discomfort. It is a reliable aid to combat sleeping disorders and other ailments, encouraging inner peace and restful sleep.

### Areas of application:

- Pain reduction
- Improvement of blood circulation
- Relaxation
- Sleeping disorders, stress
- Sore muscles

### Side effects

After the application a slight sensation of drowsiness may occasionally occur. This is completely harmless. In rare cases momentary headaches or dizziness may occur after an application. These constitute a natural reaction by the body to the effect of the SISSEL® Acupressur Mat.

Patients suffering from a whiplash may initially experience neck pain or dizziness. We recommend to initially reduce the application duration to 15 minutes every other day. Afterwards the treatments may be gradually intensified in the respective rhythm.

Patients suffering from fibromyalgia and increased sensitivity to pain should initially limit the use to a short period and use a thicker fabric layer between body and mat.

### Warning

- Do not hurt your hands when getting on or off the mat!
- Should you experience nausea or dizziness during a treatment, interrupt the application and place a thicker padding on the mat.
- People with sensitive skin should be careful and always place a towel or T-shirt on their SISSEL® Acupressur Mat in order to avoid itching.
- People suffering from disordered blood clotting should abstain from using the SISSEL® Acupressur Mat in order to avoid skin bleeding and bruising.
- People with severe allergic back eczema or general psoriasis should exercise special care and always use a towel. Initially we recommend shorter resting sessions on your SISSEL® Acupressur Mat to gradually find the optimal treatment duration.
- People with extended areas of damaged skin on their backs should not use the SISSEL® Acupressur Mat.

# ACHING MUSCLES

## Set Up

The SISSEL® Acupressur Mat is a two piece set. Attach the pieces on top of a firm surface. Total size is 40 x 60 cm. Place the SISSEL® Acupressur Mat onto a bed, a massage table or on a rug on the floor.

- The pressure can be individually regulated with the help of a towel or T-shirt.
- The first minutes on the mat might be experienced as unpleasant. This should improve significantly, however, to be followed by a sensation of well-being.
- Recommended treatment duration is about 15 - 40 minutes, depending on the needs and physical well-being (or the recommendation of your doctor or therapist.) The SISSEL® Acupressur Mat may also be used as a precaution for health promotion.
- The SISSEL® Acupressur Mat can assist you in overcoming sleeping disorders. Rest on the mat for some time until you have calmed down. Do not sleep on the mat!
- Warning! Use the mat exclusively for back treatments.

## Storage

When storing the mat avoid any unintentional contact with the spikes (eyes!).


## Cleaning

Rinse the SISSEL® Acupressur Mat with hot water and ordinary household cleaner.

## General

If you have questions about the product or the applications please talk to your physician or therapist.

- Size: 40 x 60 cm
- Adjustable effect
- Material: recyclable, hygienic TPE (can be cleaned and disinfected)

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**SISSEL® GmbH**