The natural way of Sweden

SISSEL® Spiky-Ball

- | Gebrauchsanweisung
- | Instructions for use
- | Mode d'emploi
- | Gebruiksaanwijzing
- | Istruzioni per l'uso
- | 사용방법



1.


2.


3.

4.

SISSEL® Spiky-Ball

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SISSEL® Spiky-Ball

For gentle and stimulating massage – fast, easy,
can be used anywhere

Function

The SISSEL® Spiky-Ball is a simple but effective aid to stimulate the skin which transmits a stimulus to the nerve pathways and the muscles underneath. This improves blood circulation which is noticeably felt as soothing warmth. The ball is available in different sizes and can be used in different ways. Small Spiky-Balls are more suitable in particular for gripping exercises and to massage arms, hands, and feet. Balls with a diameter of 9-10 cm are ideal for massaging the back and legs.

Range of application

The SISSEL® Spiky-Ball provides fast and effective relief from tension particularly caused by stressful situations which is often felt in the neck area and the back. a gentle massage of the affected body parts. The autonomic nervous system is harmonized and muscle tension relieved.

- massage to relieve tension and improve circulation
- self-massage, partner massage, child massage (not under the age of three)
- reflexology
- training body awareness
- gripping and strengthening exercises for fingers and hands
- exercises for varicose veins, massaging aching feet
- relaxing exercises in the office or at the computer

Instructions for use

Roll the Spiky-Ball over the body, directly onto the skin or over clothing, using small circular movements. The muscles should be relaxed. The hardness of the ball and therefore the intensity can be regulated using the valve. Massage the muscular tissue gently at first and over a large area using long rolling motions before targeting the problem areas with firmer circling movements. You can also hold the Spiky-Ball between your back and a wall to massage your back and gluteal muscles.

Safety guidelines

- The sensation of the stimulus from the Spiky-Ball massage varies in the different parts of the body (but the pressure felt should always be pleasant). Stop massaging immediately if you experience pain.
- **Do not use on bones and joints, especially not on the spine.**
- **Do not use on the eyes or on other sensitive body areas.**
- **Do not use on people with limited capacities.**
- Massage larger areas moving the ball towards the heart/center of the body.
- Consult your doctor if you have any questions.

Storage/Care

- protect from direct sunlight and intense heat
- clean using normal household cleaners, disinfectant spray or wipes

Material/Disposal

- material: phthalate-free PVC
- disposal: with household waste or use a recycling option where applicable

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SISSEL® GmbH